

Heytemp - T15

Sintering Furnace

Maximum temperature

High Purity Elema Technology

90mm*2 super-diameter tray

Up to 80 teeth can be sintered



Advantage:

- · Fast crystallization of full anatomical zirconia
- Pre-drying time can be set to completely dry the denture
- Maximum heating rate 200 degrees per minute
- It is ready to use without preheating in advance
- · Automatic fast cooling technology
- Closed kiln design, with maximum safety(chair side)
- Variable load design, automatic power regulation of sintering furnace(energy saving)



Fast temperature rise to 1530°C in 20 minutes!